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Estimated run time (minutes) including setup and treatment times. Runtimes greater than 10 hours are highlighted.

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Imported/Distributed in Canada by / Distribué/Importé au Canada par:
Fresenius Medical Care Canada Inc.
45 Staples Avenue, Suite 110, Richmond Hill, Ontario L4B 4W6
Service and complaints/Service technique et plaintes 1-888-709-4411

Imported/Distributed in Mexico by / Importado y distribuido en México por:
Fresenius Medical Care de México S.A. de C.V.
Paseo del Norte 5300, "Guadalajara Technology Park", CP 45010
Zapopan, Jal. México (52) (33) 35 40 42 00
Servicio de atención a clientes: 01. 800. 483. 8383

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**FRESENIUS
MEDICAL CARE**

RENAL THERAPIES GROUP

Powder for production of liquid sodium bicarbonate concentrate for bicarbonate hemodialysis

INDICATIONS FOR USE: The Fresenius Medical Care North America (FMCNA) *bibag* system is used with FMCNA three stream proportioning hemodialysis systems equipped with the *bibag* module, such as the 2008T Hemodialysis Machine, and is intended for use in bicarbonate hemodialysis for acute and chronic renal failure. *bibag* bicarbonate concentrate is intended for extracorporeal bicarbonate hemodialysis according to a physician's prescription.

CONTRAINDICATIONS: There are no absolute contraindications for bicarbonate dialysis. *bibag* bicarbonate concentrate must be used only on the direction of a physician who has evaluated all the pertinent features of the product in relation to the individual patient.

REQUIREMENTS:

- For use with FMCRTG dialysis machines with three-stream proportioning bicarbonate systems when set to proportion 1 part acid to 44 parts bicarbonate base and water.
- For use only with FMCRTG systems equipped with the *bibag* connector and Diasafe Plus filter.
- For use only with FMCRTG 45x acid concentrates.
- For use only with water that meets current AAMI RD62 or ISO 13959 requirements.

WARNING:

- Failure to follow these Instructions for Use may cause patient injury or death.
- Test the conductivity and approximate pH of the dialysate with an independent device before beginning treatment. Use of a dialysate with incorrect conductivity or pH can cause serious injury or death.
- The dissolved *bibag* bicarbonate concentrate must be used within 24 hours of connecting to the dialysis machine. Bacterial growth may occur when using bicarbonate concentrate.

Refer to the Hemodialysis Machine and *bibag* System Manuals for a complete list of Warnings associated with the use of the *bibag* System in hemodialysis treatment.

CAUTION:

- Federal (USA) law restricts this device to sale by or on the order of a physician. 
- Do not use if the product is open or damaged.
- Replace a leaking *bibag* immediately. Spills can cause slips and falls; clean up spills immediately.
- Product must be used before the expiration date on the product's label.

ADVERSE REACTIONS: Some undesirable effects related to dialysis can occur, for instance, hypotension, nausea, headache and cramps.

DIRECTIONS FOR USE:

Note: See table on front page for treatment times for 650g/900g *bibag* at various dialysate flow rates and sodium bicarbonate concentrations.

1. Remove the plastic seal from underneath the water and bicarbonate nozzles of the *bibag*. Then open the *bibag* door on the machine by lifting up on the dark-gray handle.
2. With the white *bibag* handle facing outward, hang the *bibag* on the *bibag* connector nozzles. Close the door, making sure it latches firmly in place.
3. Refer to the detailed directions for use provided in the hemodialysis machine and *bibag* Systems manuals.
4. At the end of treatment, use the "Empty *bibag*" function of the dialysis machine (or drain to a convenient sink) to empty the bag for easier disposal.
5. Lift up on the dark-gray *bibag* door handle to open the *bibag* door. Remove the *bibag* and dispose of it and any unused powder according to dialysis unit procedures.
 - Dispose of the bag and unused powder between each treatment.
 - In some cases the bicarbonate powder may contain particles that have been slightly discolored after exposure to water. These discolorations originate from the raw material (sodium bicarbonate). The powder from the *bibag* meets the requirements of the European and US Pharmacopoeia. The safety and efficacy of the *bibag* product are not affected by these discolorations.

STORAGE: Store between +5°C and +30°C (41°F and 86°F). Do not refrigerate. In order to avoid any contamination, keep the *bibag* tightly sealed until it is connected to the hemodialysis machine.

CHEMICAL COMPOSITION OF CONCENTRATE: Sodium Bicarbonate, USP/Ph Eur

PRESENTATION: 650g or 900g *bibag* bicarbonate concentrate in polyamide/polyethylene container.

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