
2008 Series Hemodialysis Machine Operator's Manuals Addendum for Concentrate Connection

Caution: Federal (US) law restricts this device to sale only by or on the order of a physician.

This is an addendum to the 2008 Series Hemodialysis Machine Operators Manuals.

This addendum relates to the following additions:

- Clarification on Dialysate Path and connecting the machine to acid and bicarbonate concentrate.
- Instruction on proper usage and setup of the wand assemblies with the 2008 Series Hemodialysis Machines.

Indications for Use Statement:

The Wand Assemblies are intended to be used with the 2008 Series Hemodialysis Machine to provide a fluid pathway for the concentrates and acetic acid from the canister to the acid and bicarbonate lines of 2008 Series Hemodialysis Machines. The Wand Assemblies are designed to be used as an accessory of the 2008 Series Hemodialysis Machine.

The Dialysate Path


Figure 1

Dialysate Path

The concentrate connectors draw acid and bicarbonate concentrates from the supply canister. Next the machine mixes the concentrates with AAMI Standard water into dialysate. The dialysate then runs through the dialyzer and out the machine's drain line.

Red acid and blue bicarbonate connectors

The red and blue concentrate connectors shown above in Figure 1 draw acid and bicarbonate concentrates from the supply canisters. The concentrate connectors pull out of the machine and snap into the acid and bicarbonate canister with the help of a connector cap. Without the connector cap, the acid and bicarbonate connectors would not be able to connect to the appropriate canister. The connector cap aids the transfer of acid and bicarbonate from the canister to the machine. Follow the steps on the next page in order to successfully attach the connector cap to the canister and the concentrate connector to the cap.

Concentrate Canister Connector Cap

The connector cap screws onto the top of concentrate canister. The Acid and Bicarbonate connectors snap into the caps so the machine can pull concentrate from the canisters.

Connecting the Concentrates


Figure 2

Figure 2 above demonstrates the process to properly connect the cap and the concentrate connector to the canister. The blue (bicarbonate) side shows the step by step process while the red (acid) side shows the finished product after following the proper steps.

1. Check Acid and Bicarbonate Canister

Before startup, make sure that the acid and bicarbonate canisters are at appropriate levels.

2. Attach Connector caps

Screw on the connector caps to the canister. Make sure they are tight and that the tip of the wand is submerged in concentrate.

3. Remove Concentrate Connectors

Remove the acid (red) and bicarbonate (blue) connectors from the front of the machine and attach them to the connector caps on the appropriate canister.

4. Check System

After step 3, make sure that all the connections are tight and secure. Once they are, you are ready to successfully pull acid and bicarbonate concentrates from the canister.

Cleaning

Follow your clinic's instructions on disinfecting the concentrate canister and connector caps. Refer to the 2008 Series Hemodialysis Machine Operator's Manuals for additional information.

Please call Fresenius Medical Care Technical Support with any questions at (800) 227-2572

Distributed by:
Fresenius USA, Inc.
4040 Nelson Ave
Concord, CA 94520
(800) 227-2572

THIS PAGE LEFT BLANK INTENTIONALLY

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