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CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

DIN 02325586

2360C2011 SAA

Estimated run time (minutes) including setup and treatment times.

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4D5011 / 2

00478 Rev C 11/17

Powder for production of liquid sodium bicarbonate concentrate for bicarbonate hemodialysis

INDICATIONS FOR USE: The bibag system is indicated for use in patients undergoing extracorporeal bicarbonate hemodialysis for acute and chronic renal failure. The bibag system is intended to be used as one component in the preparation of dialysate according to a physician's prescription in a 3-stream proportioning hemodialysis machine equipped with the bibag module.

CONTRAINDICATIONS: There are no absolute contraindications for bicarbonate dialysis. bibag bicarbonate concentrate must be used only on the direction of a physician who has evaluated all the pertinent features of the product in relation to the individual patient.

REQUIREMENTS:

- For use with FMCRTG dialysis machines with three-stream proportioning bicarbonate systems when set to proportion 1 part acid to 44 parts bicarbonate base and water.
- For use only with FMCRTG systems equipped with the bibag connector and Diasafe Plus filter.
- For use only with FMCRTG 45x acid concentrates.
- For use only with water that meets current AAMI RD62 or ISO 13959 requirements.

WARNING:

- Failure to follow these Instructions for Use may cause patient injury or death.
- Test the conductivity and approximate pH of the dialysate with an independent device before beginning treatment. Use of a dialysate with incorrect conductivity or pH can cause serious injury or death.
- The dissolved bibag bicarbonate concentrate must be used within 24 hours of connecting to the dialysis machine. Bacterial growth may occur when using bicarbonate concentrate.

Refer to the Hemodialysis Machine and bibag System Manuals for a complete list of Warnings associated with the use of the bibag System in hemodialysis treatment.

CAUTION:

- Do not use if the product is open or damaged. 
- Replace a leaking bibag immediately. Spills can cause slips and falls; clean up spills immediately.
- Product must be used before the expiration date on the product's label.

ADVERSE REACTIONS: Some undesirable effects related to dialysis can occur, for instance, hypotension, nausea, headache and cramps.

DIRECTIONS FOR USE:

Note: See table on front page for treatment times for 650g/900g bibag at various dialysate flow rates and sodium bicarbonate concentrations.

1. Remove the plastic seal from underneath the water and bicarbonate nozzles of the bibag. Then open the bibag door on the machine by lifting up on the dark-gray handle.
 2. With the white bibag handle facing outward, hang the bibag on the bibag connector nozzles. Close the door, making sure it latches firmly in place.
 3. Refer to the detailed directions for use provided in the hemodialysis machine and bibag Systems manuals.
 4. At the end of treatment, use the "Empty bibag" function of the dialysis machine (or drain to a convenient sink) to empty the bag for easier disposal.
 5. Lift up on the dark-gray bibag door handle to open the bibag door. Remove the bibag and dispose of it and any unused powder according to dialysis unit procedures.
- Dispose of the bag and unused powder between each treatment.
 - In some cases the bicarbonate powder may contain particles that have been slightly discolored after exposure to water. These discolorations originate from the raw material (sodium bicarbonate). The powder from the bibag meets the requirements of the European and US Pharmacopoeia. The safety and efficacy of the bibag product are not affected by these discolorations.

STORAGE: Store between +5°C and +30°C (41°F and 86°F). Do not refrigerate. In order to avoid any contamination, keep the bibag tightly sealed until it is connected to the hemodialysis machine.

CHEMICAL COMPOSITION OF CONCENTRATE: Sodium Bicarbonate, USP/Ph Eur.

PRESENTATION: 650g or 900g bibag bicarbonate concentrate in polyamide/polyethylene container.

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